

EXPERIMENTAL WORK ON WATER IN TOMSK

L. SINITSIA

Institute of Atmospheric Optics SB RAS,
1, Akademicheskii Av., Tomsk, 634055, Russia

Experimental investigations of water vapor absorption spectra at the Institute of Atmospheric Optics are performed using high sensitive Intracavity Laser Spectrometers and Grating Spectrophotometers in high frequency region $9000 - 14000 \text{ cm}^{-1}$. This region is occupied by transitions to high resonance polyads involving three or more vibrational quanta. Corresponding energy levels lie near or above barrier to linearity – saddle point on the potential energy surface, which determines in a great extend the vibration – rotation dynamics in the XH_2 – type molecules.

The high temperature water vapour absorption spectra measurements

Measurements of the high temperature water vapour absorption spectra (up to $T=800 \text{ K}$) in the region of $10000 - 16000 \text{ cm}^{-1}$ will be made using grating spectrometer.

Parameters of the spectrometer:

Focal length	2600 cm
Spectral resolution	0.03 cm^{-1}
Absorption pathlength	1-2 m

Measurements of line position and strengths, lines broadening and shifting coefficients are producing at the spectral resolution 0.03 cm^{-1} and at different sample pressures. The main advantage is the use of the diode array detector of high sensitivity that allows reaching the high signal – to – noise ratio from 1000 to 5000 depending on wavelength.

Intracavity laser spectrometers

The H_2O absorption spectra measurements are performed using Intracavity laser spectrometers operating in the region upper 8500 cm^{-1} .

Lasers for ICL-spectroscopy:

1. Nd-glass ($9200-9520 \text{ cm}^{-1}$)
2. Color center in LiF ($10600-11200 \text{ cm}^{-1}$)
3. Ti:Sapphire ($11300-12800 \text{ cm}^{-1}$)

Lasers have pulse duration 25-2000 microseconds that allow getting high threshold sensitivity to absorption.

The measurements will be made at spectral resolution 0.03 cm^{-1} and effective path length up to 10 km. Spectra of the H_2^{16}O , HDO , D_2O and H_2^{18}O isotopomers will be recorded at elevated temperatures, up to 800 K and sample pressures up to 100 Torr.